Standard Charger

MEENERGY™ BTE P/SP Charger Manual

www.miracle-ear.com

Content

Before you start **4**

Components 6

Connecting the USB power cable 7

Charging **8**

Correct position of the hearing aids 10

Charging status 11

Charging tips 12

Further information **19**

Maintenance and care 20

Troubleshooting 20

Important safety information **23**

General warnings 23

Personal safety 25

Product safety 30

Important information **32**

Explanation of symbols 32

Operating, transport, and storage conditions 34

Technical data of the power supply adapter 36

Disposal information 36

Conformance information 37

Before you start

This device is intended to charge our specific hearing aids with built-in power cells (lithium-ion rechargeable batteries). Your Hearing Care Specialist will advise you on the compatible models.

The charger is only intended for indoor use. Observe and follow the operating conditions in section "Operating, transport, and storage conditions".


The charger is not intended as a means of storage for transportation of hearing aids.

CAUTION

Read this user guide thoroughly and completely and follow the safety information in this document to avoid damage or injury.


Components

-
- | | |
|------------------|----------------------------------|
| 1 Status LEDs | 3 Non-detachable USB power cable |
| 2 Charging slots | 4 Power adapter |
-

Connecting the USB power cable

- ▶ Connect your USB power cable to the provided power adapter.

Charging

NOTICE

Please make sure that the power plug is easily accessible to remove it from power supply if necessary.

- ▶ Connect the charger to the power supply. LEDs for each charging slots blink 3 times to indicate that the charger has power and is ready for use.

- ▶ Put the hearing aids in the charging slots as shown in the picture. You can charge both hearing aids at the same time. The hearing aids are automatically turned off and start charging. The charging status is indicated, see table "Charging status" on the next pages.

- ▶ When you want to use your hearing aids, take them out of the charger. They will automatically turn on.
- ▶ You can leave the charger plugged in or remove it from the power supply.

Correct position of the hearing aids

- ▶ If not indicated differently, the left and right hearing aids can be inserted into either charging slot. The slots are not side specific.
- ▶ Make sure the bottom parts of the hearing aids are placed well into the charging slots. Correct placement will be indicated by the LED.

Charging status

LED	Status description
	No LED light: Charger is not connected to power supply. No hearing instrument is inserted in the respective charging slot or the hearing instrument is not inserted correctly.
	Green LED light blinks 3 times: Charger is connected to a power supply and ready to charge.
	Flashing green LED light: Charging in progress.

LED	Status description
	Solid green LED light: Hearing instrument is fully charged.
	Red LED light blinks 3 times in 2 seconds: Error Read more in section "Troubleshooting".

Charging tips

You may have some questions on how to best charge your hearing aids. Read the questions and answers below. If you have more questions, your Hearing Care Specialist will help you.

How often do I need to charge the hearing aids?

We recommend charging the hearing aids daily, even if the power cell is not empty. Charge them overnight so you can be sure to start your day with fully-charged hearing aids.

How long does charging take?

A full charging cycle takes approximately 3 to 4 hours. After 30 minutes, it is possible to use your hearing aids for up to 7 hours.

Charging times may vary due to the use of the hearing aids and the age of the power cell.

Can I leave the hearing aids in the charger?

Even when fully charged, you can leave the hearing aids in the charger. They automatically stop charging when they are fully charged.

What happens if I forgot to charge?

Your hearing aids will automatically switch-off when the power cell is empty. Charge the hearing aids as soon as possible.

Hearing aids which have not been charged for over 6 months may be damaged due to an unrecoverable deep-discharge of the power cell. Deep-discharged power cells cannot be charged anymore and have to be replaced.

What if I intend not to use the hearing aids for several days or longer?

- Short-term storage (up to several days): Turn off your hearing aids via the push button or place them in the powered charger.

The charger must be connected to the power supply. If the charger is not connected to power supply and you place your hearing aids in the charger, the hearing aids are not turned off.

Note that when you turn off the hearing aids via a remote control or smartphone app, the hearing aids are not turned off completely. They are in standby-mode and still draw some power.

- Long-term storage (weeks, months, ...): First, fully charge your hearing aids. Turn them off via the push button before storing them. We recommend to use drying products while storing the hearing aids. At 6 month intervals you will have to charge the hearing aids to avoid an unrecoverable deep-discharge of the power cell. Deep-discharged power cells cannot be charged anymore and they have to be replaced. We recommend recharging more frequently than 6 months intervals. Observe the storage conditions in hearing aids' user guide. Ask your Hearing Care Specialist for recommended drying products.

What should I pay attention to when charging?

The hearing aids and charger must be clean and dry. Observe the operation conditions in section "Operating, transport, and storage conditions". Avoid charging at high ambient temperatures, as this can shorten the life span and capacity of the power cell.

If you use the charger in a very warm environment, note that it will stop charging in order to cool down if the charger's temperature exceeds 113° F. The LED light will blink 3 times in 2 seconds when this happens. Charging resumes automatically after the device has cooled down. Do not touch the hearing aids, they may be hot.

What should I pay attention to when using a power strip with energy-saving features?

- Turn on the power strip before you insert the hearing aids in the charger. Otherwise, the hearing aids may not be detected by the charger.
- Do not turn off the power strip while the hearing aids are in the charger. Otherwise the hearing aids will automatically turn on and discharge in the charger. Before you turn off the power strip, remove the hearing aids from the charger and turn them off via the push button.

Further information

Maintenance and care

NOTICE

Risk of injuries caused by electric shock!

- ▶ If cleaning is required, disconnect the charger from the power supply.
- ▶ Clean the charger with a soft, dry tissue.
- ▶ Never use running water or immerse the charger in water.
- ▶ Do not sterilize or disinfect the charger.

- ▶ After use and before placing in the charger, wipe the charging pins on the hearing aids with a soft, dry tissue to remove grease or moisture.

Troubleshooting

Problem and possible solutions

LED light is off.

- Check if the power supply is connected correctly.

LED light is blinking red 3 times.

An error in the charger has occurred.

- Take the hearing aids out of the charger and disconnect the charger from the power supply.
 - Reconnect the charger to the power supply and re-insert the hearing aids in the charging slots.
 - Charger has reached the temperature limit and needs to cool down. Charging is resumed when the temperature decreases below 113° F.
-

If the LEDs are still blinking red:

- Take the hearing aids out of the charger and disconnect the charger from the power supply.
- Contact your Hearing Care Specialist for support.

Consult your Hearing Care Specialist if you encounter further problems.

Important safety information

General warnings

- Observe the operating conditions. Refer to the section “Operating, transport, and storage conditions”.
- Pay special attention to protect your device from extreme heat (e.g. fire, microwave oven, induction oven or other high induction fields) and sunlight.
- Charge your hearing aids only with the approved charger. Ask your Hearing Care Specialist for support.
- Do not disassemble the device.
- Do not use the device in explosive atmospheres.
- Do not use devices that are deformed or obviously damaged.

- Do not use the device when its performance after charging decreases considerably. Return the device to your Hearing Care Specialist.
- Do not dispose of the device in fire or water.

In any of the above cases, return your device to the Hearing Care Specialist.

WARNING

Risk of electric shock!

- ▶ Connect the charger only with a USB power supply which complies with IEC 60950-1 and/or IEC62368-1.

Personal safety

WARNING

Risk of burns!

- ▶ Do not place anything other than your hearing aids in the charger.

WARNING

Risk of electrical shock!

- ▶ Do not use obviously damaged devices and return them to point of sale.

WARNING

Choking hazard!

Your device contains small parts which can be swallowed.

- ▶ Keep hearing aids, batteries and accessories out of reach of children and mentally disabled persons.
- ▶ If parts have been swallowed consult a physician or hospital immediately.

WARNING

Note that any unauthorized changes to the product may cause damage to the product or cause injury.

- ▶ Use only approved parts and accessories. Ask your Hearing Care Specialist for support.

WARNING

Risk of explosion!

- ▶ Do not use your device in explosive atmospheres (e. g. in mining areas).

WARNING

Charging systems may interfere with measuring devices and electronic equipment.

- ▶ In hospitals: The charger complies with CISPR 35 and CISPR 32. Therefore the use of the charger in hospitals should be permissible. Nevertheless the charger may interfere with highly sensitive devices. In those cases appropriate steps shall be taken (increase distance from sensitive devices or disconnect the charger from power supply).

WARNING

Interference with implantable devices

- ▶ The device can be used with all electronic implants that comply with the ANSI/AAMI/ISO 14117:2012 Electromagnetic compatibility standard for human implants.
- ▶ In case your implant does not comply with the ANSI/AAMI/ISO 14117:2012 standard, consult the manufacturer of your implantable device for information about the risk of disturbance.

Product safety

NOTICE

- ▶ Protect your devices from extreme heat. Do not expose them to direct sunlight.

NOTICE

- ▶ Protect your devices from high humidity.

NOTICE

- ▶ Do not dry your devices in the microwave oven.

NOTICE

- ▶ Do not drop. Internal components may get damaged.

NOTICE

- ▶ Do not disassemble the device.

NOTICE

- Use of incompatible parts poses a risk of malfunction.
- ▶ Use only the power supply and country specific adapter that is provided with the charger.

Important information

Explanation of symbols

Symbols used in this document

-  Points out a situation that could lead to serious, moderate, or minor injuries.
-  Indicates possible property damage.
-  Advice and tips on how to handle your device better.

Symbols on the device or packaging

-  CE compliance label, confirms compliance with certain European Directives, refer to section "Conformance information".
-  Indicates the legal manufacturer of the device.
-  Do not dispose of the device with general domestic waste. Read more in section "Disposal information".
-  Read and follow the instructions in the user guide.

Symbols on the device or packaging

Medical device

Operating, transport, and storage conditions

Operating conditions

Temperature	50 to 95° F
Relative humidity	20 to 75%

Operating conditions

Atmospheric pressure 700 to 1060 hPa

During extended periods of transport and storage, please observe the following conditions:

	Storage	Transport
Temperature	50 to 104° F	-4 to 140° F
Relative humidity	10 to 80%	5 to 90%
Atmospheric pressure	700 to 1060 hPa	700 to 1060 hPa

For other parts, such as hearing aids, other conditions may apply.

Technical data of the power supply adapter

Input voltage:	100-240 V AC
Input current:	0.5 A
Input frequency:	50-60 Hz
Output voltage:	5 V DC
Output current:	1 A

Disposal information

Recycle devices, accessories and packaging according to local regulations.

To avoid environmental pollution, do not throw the device into household trash. Dispose of devices according to local regulations.

Conformance information

This device complies with Part 15B of the FCC Rules and with CAN ICES-003 / NMB-003.

Operation is subject to the following conditions:

- This device may not cause harmful interference,
- this device must accept any interference received, including interference that may cause undesired operation, and
- the user shall keep a minimal distance of 20 cm (8 in) during the charging process.

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Copyright ©2021 Miracle-Ear, Inc.
Distributed By Miracle-Ear, Inc.
Derechos de autor ©2021 Miracle-Ear, Inc.